



Model	Maximum operating temperature (°C)	Capacity (l)	Internal dimensions (mm – h x Ø)	External dimensions (mm – h x w x d)	Heat up time to 100°C below maximum temperature (minutes)	Maximum operating power (kW)
EHF 17/3	1700	3.4	190 X 150	975 X 750 X 530	80	8
EHF 17/8	1700	7.9	250 X 200	1950 X 1360 X 800	80	9
EHF 18/3	1800	3.4	190 X 150	975 X 750 X 530	110	6
EHF 18/8	1800	7.9	250 X 200	1950 X 1360 X 800	110	9