



Model	400°C 500°C 600°C	WHT 4/30 WHT 5/30 WHT 6/30	WHT 4/60 WHT 5/60 WHT 6/60	WHT 4/120 WHT 5/120 WHT 6/120
Chamber dimensions (mm – h x w x d)		300 x x 300 x 305	400 x 400 x 405	650 x 480 x 405
External dimensions (mm – h x w x d)		570 x 830 x 570	670 x 930 x 670	920 x 1030 x 670
Usable volume (litres)		30	60	120
Heat up time @ 240V (minutes)	400°C 500°C 600°C	50 75 120	50 75 120	50 75 120
Recovery time (minutes) 240V	400°C 500°C 600°C	10 16 20	10 16 20	10 16 20
Temperature stability: PID control		± 0.5°C	± 0.5°C	± 0.5°C
Temperature uniformity measured @ 600°C		± 5°C	± 5°C	± 5°C
Power (kW)	400°C 500°C 600°C	1.0 2.0 2.0	1.5 2.25 2.25	2.25 3.0 3.0
Number of shelves supplied		2	2	2
Weight (kg)		73	99	179