



Model	Maximum operating temperature (°C)	Maximum continuous temperature (°C)	Tube inside diameter (mm)	Heated length (mm)	Maximum power (kW)	External dimensions including console (mm – h x l x d)	Optional separate work tubes (mm – id x length)	
							Short	Long
CSC 12/90/300 (H)	1200	1150	75	300	1.5	350 x 425 x 410	75 x 450	75 x 750
CSC 12/90/300 (V)	1200	1150	75	300	1.5	650 x 350 x 360	75 x 450	75 x 750
CSC 12/90/450 (H)	1200	1150	75	450	2.25	350 x 575 x 410	75 x 600	75 x 900
CSC 12/90/450 (V)	1200	1150	75	450	2.25	800 x 350 x 360	75 x 600	75 x 900
CSC 12/90/600 (H)	1200	1150	75	600	3.0	350 x 725 x 410	75 x 750	75 x 1050
CSC 12/90/600 (V)	1200	1150	75	600	3.0	950 x 350 x 360	75 x 750	75 x 1050
CSC 17/90/250 (V)	1700	1650	75	250	4.5	865 x 600 x 705	75 x 550	75 x 900